

DW1331 SIEVE ANALYSIS TEST APPARATUS



Model No.

DW1331

Application:

DW1331 sieve analysis test apparatus is used to assess the particle size distribution (also called gradation) of granular or powder material. It is suitable for particle gradation test, impurity content test and liquid filtration test.

Related Standards:

STANDARDS

Note: this tester can conform to but not limit to all the standards above, for more standards conformance, please contact us.

Features:

- Low noise, light weight, high efficiency, high precision
- Comprehensive function, easy to operate
- Particle can be accurately divided into 1~8 grades in one test
- New special power generator, high frequency vibration
- Standard sieves made in military industrial enterprises, high quality



- Three display modes of timer (seconds, minute and hour) are selectable



Key Specifications:

- | | |
|-----------------------------------|---|
| ● Inner diameter of sieves | Φ200/100/75mm |
| ● Diameter of sieve pore (Φ200) | 63μm, 90μm, 125μm, 180μm, 250μm,
355μm, 500μm, 710μm, 1000μm |
| ● Number of sieves | 8pcs |
| ● Superimposed height of sieves | 400mm |
| ● Amplitude | ≤5mm |
| ● Timer range | 0~60min |
| ● Diameter of particle | 0.038mm~3mm |
| ● Noise | ≤50dB |
| ● Feeding quantity (one time) | ≤200g |
| ● Power supply | 220V, 50Hz 125W |
| ● Motor speed | 1420 r/min |
| ● Dimension | 360×300×736mm (L×W×H) |
| ● Weight | 25kg |

Standard Configuration:

No.	Item	Quantity
-----	------	----------



WWW.FYICHINA.COM.CN

FANYUAN INSTRUMENT (HF) CO., LTD.
DOWELL SCIENCE AND TECHNOLOGY (HK) CO., LTD.
ADD: Tianyuan Road 9#, High-tech Zone, Hefei City, China Postal Code: 230088
TEL: +86-551-65232508
FAX: +86-551-65232507
Email: sales@hkcdowell.com
Website: http://www.fyichina.com.cn

1	Main Machine	1set
2	Sieves	8pcs
3	Power Cord	1pc

Optional Accessories:

No.	Item	Quantity